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The screenshot shows the 'Sieving' software interface with the 'Sample' tab selected. The interface includes a search bar, navigation tabs (Fast sieving, Manual sieving, Methods), and a sidebar with options like Basic, Sample, Methodology, Data, and Device. The main area is titled 'Define sample information' and contains the following fields:

- Characteristics:**
 - Sample material*: APHAZ
 - Sample preparation*: No sample preparation
 - Density: 100 g/cm³
 - Source: Container
 - Sampling: [empty]
 - Sample weight: 5000 g
- Sample tolerance and check:**
 - On ☒ Weigh-in tolerance: 0.00 % To 5.00 %
 - Off ☐ Weighing assistant
- Comment:** [empty text box]

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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The screenshot shows the 'Sieving' software interface with the 'Sample' tab selected. The interface is organized into several sections: 'Sample preparation' with fields for 'Dividing by' (set to 'Divider'), 'Density' (100 g/cm³), 'Source' (set to 'Container'), 'Sampling' (set to 'With Spoon'), and 'Sample weight' (10.00 g). Below this is the 'Sample tolerance and check' section, which has two options: 'On Weigh-in tolerance' (selected) and 'Off Weighing assistant'. The 'Backweighing tolerance' section is highlighted in green and contains two entries: 'Mean size (0.02 mm)' with 'Expected overgrain' set to 7.50 g and 'From'/'To' percentages of -15.00 % to 15.00 %, and another entry with 'Expected overgrain' set to 4.00 g and 'From'/'To' percentages of -15.00 % to 15.00 %. A 'Comment' field is located at the bottom of the interface.

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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TREND ANALYSIS: TEST EQUIPMENT AND PROCESS MONITORING - PROACTIVELY AVOID INCORRECT SIEVING

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
Trend Analysis with downward outliers of 5 consecutive sieve analyses (blinding)

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Search

+ Add method

Name	Sample material	Actions
Coffee Arabica	Ground Coffee	
Coffee Barista	Ground Coffee	
LAST USED	Ground Coffee	

1/11 OPTIONAL SETUP
CHOOSE A METHOD FOR EDITING

Sieving

Fast sievingManual sievingMethods

SearchCancel

BasicSampleMethodologyDataDevice

Define sample information

Characteristics

Sample materialGround Coffee

Sample preparationNo sample preparation

Density0.37 g/cm³

SourceCoffee Mill 04

Sampling

Sample weight0.00 g

Sample tolerances and check

OffWeigh-in tolerance-5.00 %To5.00 %

OnWeighing assistant

Comment

**2/11 OPTIONAL SETUP
ACTIVATION OF THE WEIGHING ASSISTANT**

Sieving

Fast sieving

Manual sieving

Methods

Search

+ Add method

Name	Sample material	Actions
Coffee Arabica	Ground Coffee	
Coffee Barista	Ground Coffee	
LAST USED	Ground Coffee	

3/11 STARTING THE METHOD

Step 1 / 7 Sieve 1 / 7

Specify the
sample

Method name

Coffee Arabica

Sample name

Arabica

Sample ID

06

Comment

Username

GBE

Department

Ok

4/11 GUIDED SIEVING
STEP 1: SAMPLE IDENTIFICATION

Step 2 / 7 Sieve 1 / 7

Remove sieve
and lid

Ok

5/11 GUIDED SIEVING
STEP 2: PREPARATION OF BALANCE

Step 3 / 7 Sieve 1 / 7

Insert the
sieve - 125 µm

205.58 g

Please provide a plausible sieve weight
>=100 g

Ok

6/11 GUIDED SIEVING
STEP 3: WEIGHING THE REQUESTED SIEVE

Step 4 / 7 Sieve 1 / 7

10.33 g

Please provide a plausible sample weight

$\geq 0.1 \text{ g}$

Weighing assistant

4.63 g - 13.62 g

Ok

7/11 GUIDED SIEVING
STEP 4: WEIGHING THE SAMPLE. THE WEIGHING ASSISTANT RECOMMENDS A SUITABLE SAMPLE AMOUNT

Step 5 / 7 Sieve 1 / 7

Put on the lid

326.13 g

Please provide a plausible lid weight

≥ 200 g

Ok

8/11 GUIDED SIEVING
STEP 5: WEIGHING THE LID

Step 6 / 7 Sieve 1 / 7

Ready for
sieving

Duration (mm:ss)

03:00

Pressure (Pa)

2,500 ^
v

Speed of nozzle drive (rpm)

20.00 ^
v

Start sieving

9/11 GUIDED SIEVING
STEP 6: OPTIONAL ADJUSTMENT OF THE SIEVING PARAMETERS AND START

Step 7 / 7 Sieve 1 / 7

Remove sieve
and lid

Please provide a plausible lift-off weight

$\leq -531.65 \text{ g}$

Next

10/11 GUIDED SIEVING
STEP 7: BACKWEIGHING THE RESIDUAL MASS

Preliminary results

125 µm

Sieving completed successfully

Results

Oversize	7.73 g
----------	--------

Sample information

Sample amount	10.32 g
---------------	---------

Undersize	2.59 g
-----------	--------

Device parameters

Duration	01:30 mm:ss
----------	-------------

Speed of nozzle drive	20 rpm
-----------------------	--------

Pressure	3,000 Pa
----------	----------

Average pressure	2,784 Pa
------------------	----------

Pressure std. dev. sigma	877 Pa
--------------------------	--------

Pressure variance	768,927 Pa²
-------------------	-------------

Next

11/11 DISPLAY OF RESULTS

AIR JET SIEVING MACHINES AS 200 JET PRO / AS 200 JET PHARMA

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Graphical User Interface, guided sieving	準拠
Mode of operation	Fast sieving, manual sieving and methods
Languages	German and English
Parameter units changeable	Size: mm/µm/mesh/inch Mass: g/kg Time: s/min Pressure: mbar/Pa/psi Volume: ml/cl/l/cm3
SOPs	Yes, integrated in methods, more than 40
Individual sample labelling	Yes, for each sample, reference to method during data choice and in overview
Parameter setting per test sieve	準拠
LIMS ready	準拠
possible export format	.xls, .xlsx, .xps, .csv, .pdf, .png, .jpg, .txt, .rtf, html, .mht
Teamviewer online support	準拠
Screenshot function	準拠
Back-up function	準拠
Automatic storage of results	Yes, in machine, auto report function, PDF document or to printer

www.retsched.com/AS200jet-pro-pharma

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AIR JET SIEVING MACHINES AS 200 JET PRO / JET PHARMA

Air Jet Sieving Machines for test sieves 203 mm (8") Ø, with automatic vacuum regulation and internal balance, incl. lids for test sieves 1" and 2" height and rubber mallet - please order test sieve and MANDATORY vacuum cleaner (Item No. 22.748.0020 or Item No. 22.748.0021) separately!

30.034.0001		AS 200 100–240 V, 50/60 Hz jet pro	incl. test report acc. to EN 10204 2.2 evaluation software installed (no EasySieve required)
30.034.0002		AS 200 100–240 V, 50/60 Hz jet pharma (no EasySieve required)	incl. test report acc. to EN 10204 2.2 evaluation software installed incl. AuditTrail & User Management

VACUUM CLEANER AS 200 JET PRO / JET PHARMA

22.748.0020	Industrial vacuum cleaner ATTIX 33 230V 50/60Hz (other electrical versions available) MANDATORY vacuum cleaner for 30.034.0001 & 30.034.0002 NOT SUITABLE for 32.935.0008
32.524.0017	HEPA filter for ATTIX 33
32.524.0018	Cover for HEPA filter ATTIX 33
32.524.0019	Replacement filter bags for ATTIX 33, 5 pcs.

CYCLONE AS 200 JET PRO / JET PHARMA

22.935.0042	High-performance cyclone including holder, collecting receptacle (250 ml and 500 ml sample bottle) and 1.5 m FDA certified hose with clamps (industrial vacuum cleaner 22.748.0020 or 22.748.0021 recommended)
05.707.0208	FDA certified hose, 1.5 m
02.183.0151	Connector for cyclone

ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0047

IQ/OQ Documentation for AS 200 jet pro

99.200.0046

IQ/OQ Documentation for AS 200 jet pharma

RECALIBRATION AIR JET SIEVING MACHINES

The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents

SER-97.001.0021

Calibration of sieving machines AS 200 jet pro & AS 200 jet pharma
at Retsch/Germany, incl. test report acc. to EN 10204 2.2