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The screenshot shows the 'Sieving' software interface with the 'Sample' tab selected. The interface includes a search bar, navigation tabs (Fast sieving, Manual sieving, Methods), and a sidebar with options like Basic, Sample, Methodology, Data, and Device. The main area is titled 'Define sample information' and contains the following fields:

- Characteristics:**
  - Sample material\*: APHAZ
  - Sample preparation\*: No sample preparation
  - Density: 100 g/cm³
  - Source: Container
  - Sampling: [empty]
  - Sample weight: 5000 g
- Sample tolerance and check:**
  - On ☒ Weight-in tolerance: -5.00 % To 5.00 %
  - Off ☐ Weighing assistant
- Comment:** [empty text box]

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The screenshot shows the 'Sieving' software interface with the 'Sample' tab selected. The interface includes a search bar, navigation buttons (Back, Forward, Home, etc.), and a sidebar with tabs: Basic, Sample, Methodology, Data, and Device. The 'Sample' tab contains the following settings:

- Sample preparation:** Dividing by (Dropdown), Density (100 g/cm³), Source (Container), Sampling (With Spoon), Sample weight (10.00 g).
- Sample tolerance and check:** On (Weigh-in tolerance) with a range from -5.00 % to 5.00 %, and Off (Weighing assistant).
- Screening tolerance:** On (Mean size [0.02 mm]) with an expected overgrain range from -15.00 % to 15.00 %.
- Mean size [0.02 mm]:** On (Expected overgrain) with a range from -15.00 % to 15.00 %.

At the bottom, there is a 'Comment' field.

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
Search

+ Add method

| Name           | Sample material | Actions                                                             |
|----------------|-----------------|---------------------------------------------------------------------|
| Coffee Arabica | Ground Coffee   | <div><div>▶</div><div><div></div><div></div><div></div></div></div> |
| Coffee Barista | Ground Coffee   | <div><div>▶</div><div><div></div><div></div><div></div></div></div> |
| LAST USED      | Ground Coffee   | <div><div>▶</div><div><div></div><div></div><div></div></div></div> |

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Characteristics

Sample materialGround Coffee

Sample preparationNo sample preparation

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Search

+ Add method

| Name           | Sample material | Actions                           |
|----------------|-----------------|-----------------------------------|
| Coffee Arabica | Ground Coffee   | <div><div></div><div></div></div> |
| Coffee Barista | Ground Coffee   | <div><div></div><div></div></div> |
| LAST USED      | Ground Coffee   | <div><div></div><div></div></div> |

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Insert the  
sieve - 125 µm

205.58 g

*Please provide a plausible sieve weight*  
*>=100 g*

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Step 4 / 7 Sieve 1 / 7

10.33 g

Please provide a plausible sample weight

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Step 7 / 7 Sieve 1 / 7

Remove sieve  
and lid

*Please provide a plausible lift-off weight*

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### SOFTWARE FEATURES

|                                          |                                                                                                                |
|------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| Graphical User Interface, guided sieving | 是                                                                                                              |
| Mode of operation                        | Fast sieving, manual sieving and methods                                                                       |
| Languages                                | German and English                                                                                             |
| Parameter units changeable               | Size: mm/µm/mesh/inch<br>Mass: g/kg<br>Time: s/min<br>Pressure: mbar/Pa/psi<br>Volume: ml/cl/l/cm <sup>3</sup> |
| SOPs                                     | Yes, integrated in methods, more than 40                                                                       |
| Individual sample labelling              | Yes, for each sample, reference to method during data choice and in overview                                   |
| Parameter setting per test sieve         | 是                                                                                                              |
| LIMS ready                               | 是                                                                                                              |
| possible export format                   | .xls, .xlsx, .xps, .csv, .pdf, .png, .jpg, .txt, .rtf, html, .mht                                              |
| Teamviewer online support                | 是                                                                                                              |
| Screenshot function                      | 是                                                                                                              |
| Back-up function                         | 是                                                                                                              |
| Automatic storage of results             | Yes, in machine, auto report function, PDF document or to printer                                              |

[www.retsch.cn/AS200jet-pro-pharma](http://www.retsch.cn/AS200jet-pro-pharma)

## ORDER DATA

### AIR JET SIEVING MACHINES AS 200 JET PRO / JET PHARMA

**Air Jet Sieving Machines for test sieves 203 mm (8") Ø, with automatic vacuum regulation and internal balance, incl. lids for test sieves 1" and 2" height and rubber mallet - please order test sieve and MANDATORY vacuum cleaner (Item No. 22.748.0020 or Item No. 22.748.0021) separately!**

|             |                                                                                   |                                                                        |                                                                                                               |
|-------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| 30.034.0001 |  | AS 200 100–240 V, 50/60 Hz<br>jet pro                                  | incl. test report acc. to EN 10204 2.2 evaluation<br>software installed (no EasySieve required)               |
| 30.034.0002 |  | AS 200 100–240 V, 50/60 Hz<br>jet<br>pharma<br>(no EasySieve required) | incl. test report acc. to EN 10204 2.2<br>evaluation software installed incl. AuditTrail &<br>User Management |

### VACUUM CLEANER AS 200 JET PRO / JET PHARMA

|             |                                                                                                                                                                              |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 22.748.0020 | Industrial vacuum cleaner ATTIX 33 230V 50/60Hz (other electrical versions available) MANDATORY vacuum cleaner for 30.034.0001 & 30.034.0002<br>NOT SUITABLE for 32.935.0008 |
| 32.524.0017 | HEPA filter for ATTIX 33                                                                                                                                                     |
| 32.524.0018 | Cover for HEPA filter ATTIX 33                                                                                                                                               |
| 32.524.0019 | Replacement filter bags for ATTIX 33, 5 pcs.                                                                                                                                 |

### CYCLONE AS 200 JET PRO / JET PHARMA

|             |                                                                                                                                                                                                                   |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 22.935.0042 | High-performance cyclone including holder, collecting receptacle (250 ml and 500 ml sample bottle) and 1.5 m FDA certified hose with clamps<br>(industrial vacuum cleaner 22.748.0020 or 22.748.0021 recommended) |
| 05.707.0208 | FDA certified hose, 1.5 m                                                                                                                                                                                         |
| 02.183.0151 | Connector for cyclone                                                                                                                                                                                             |

### ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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22.523.0001

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22.523.0002

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## TEST SIEVES WITH ELECTROFORMED SHEET

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## ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0047

IQ/OQ Documentation for AS 200 jet pro

99.200.0046

IQ/OQ Documentation for AS 200 jet pharma

## RECALIBRATION AIR JET SIEVING MACHINES

**The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents**

SER-97.001.0021

Calibration of sieving machines AS 200 jet pro & AS 200 jet pharma  
at Retsch/Germany, incl. test report acc. to EN 10204 2.2