

Trade name: Tungsten**Current version :** 4.0.2, issued: 03.04.2023**Replaced version:** 4.0.1, issued: 27.10.2022**Region:** GB**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1 Product identifier****Trade name****Tungsten**

Substance name Tungsten

Identification numbers

CAS no. 7440-33-7

EC no. 231-143-9

1.2 Relevant identified uses of the substance or mixture and uses advised against**Relevant identified uses of the substance or mixture**

chemical analysis

Uses advised against

No data available.

1.3 Details of the supplier of the safety data sheet**Address**

ELTRA GmbH

Retsch-Allee 1-5

42781 Haan

GERMANY

Telephone no. +49 (0) 2104 23 33-400

Fax no. +49 (0) 2104 23 33-499

Information provided by / telephone

+49 (0) 2104 23 33-192

Advice on Safety Data Sheet

info@eltra.com

1.4 Emergency telephone number

For medical advice (in German and English):

+49 (0)551 192 40 (Giftinformationszentrum Nord)

SECTION 2: Hazards identification**2.1 Classification of the substance or mixture****Classification information**

This product does not meet the classification and labelling criteria given in the Regulation (EC) No 1272/2008 (CLP).

2.2 Label elements

Not relevant

2.3 Other hazards

PBT assessment

The study does not need to be conducted according to Annex XIII of Regulation (EC) 1907/2006 (REACH).

vPvB assessment

The study does not need to be conducted according to Annex XIII of Regulation (EC) 1907/2006 (REACH).

SECTION 3: Composition/information on ingredients**3.1 Substances****Chemical characterization**

Substance name Tungsten

Formula W

Molecular weight 60.08

Identification numbers

Trade name: Tungsten**Current version :** 4.0.2, issued: 03.04.2023**Replaced version:** 4.0.1, issued: 27.10.2022**Region:** GB

CAS no. 7440-33-7
EC no. 231-143-9

3.2 Mixtures

Not applicable. The product is not a mixture.

3.3 Other information

This product is an article in the sense of article 3 of the European Chemicals Regulation (REACH) 1907/2006/EC and therefore not subject to the labelling obligations according to the EU regulation 1272/2008 (CLP). Submission of a safety data sheet is not mandatory for articles and is done on a voluntary basis.

SECTION 4: First aid measures**4.1 Description of first aid measures****General information**

In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes and launder thoroughly before reusing.

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air.

After skin contact

When in contact with the skin, clean with soap and water.

After eye contact

Remove contact lenses. Rinse eye thoroughly under running water keeping eyelids wide open and protecting the unaffected eye (at least 10 to 15 minutes).

After ingestion

Rinse the mouth thoroughly with water. Do not induce vomiting. Never give anything by mouth to an unconscious person.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically

SECTION 5: Firefighting measures**5.1 Extinguishing media****Suitable extinguishing media**

Metal fire powders; Cement; Extinguishing measures to suit surroundings.

Unsuitable extinguishing media

Water; Carbon dioxide

5.2 Special hazards arising from the substance or mixture

Combustible.

5.3 Advice for firefighters

Use self-contained breathing apparatus. Wear protective clothing.

SECTION 6: Accidental release measures**6.1 Personal precautions, protective equipment and emergency procedures****For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Ensure adequate ventilation.

For emergency responders

Personal protective equipment (PPE) - see section 8.

6.2 Environmental precautions

Do not allow to enter drains or waterways.

6.3 Methods and material for containment and cleaning up

Collect mechanically. Send in suitable containers for recovery or disposal.

Trade name: Tungsten**Current version :** 4.0.2, issued: 03.04.2023**Replaced version:** 4.0.1, issued: 27.10.2022**Region:** GB**6.4 Reference to other sections**

Information regarding safe handling, see section 7. Information regarding personal protective measures, see section 8.
Information regarding waste disposal, see section 13.

SECTION 7: Handling and storage**7.1 Precautions for safe handling****Advice on safe handling**

Provide exhaust ventilation if dust is formed. No special measures necessary if used correctly.

General protective and hygiene measures

Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages. Wash hands before breaks and after work.

Advice on protection against fire and explosion

No special measures necessary.

7.2 Conditions for safe storage, including any incompatibilities**Technical measures and storage conditions**

Keep container tightly closed and dry in a cool, well-ventilated place.

Requirements for storage rooms and vessels

No special measures required.

Incompatible products

None known

7.3 Specific end use(s)

No data available.

SECTION 8: Exposure controls/personal protection**8.1 Control parameters****Occupational exposure limit values**

No	Substance name	CAS no.	EC no.
1	Tungsten	7440-33-7	231-143-9
List of approved workplace exposure limits (WELs) / EH40			
Tungsten & compounds (as W) soluble compounds			
	WEL short-term (15 min reference period)	3	mg/m ³
	WEL long-term (8-hr TWA reference period)	1	mg/m ³
List of approved workplace exposure limits (WELs) / EH40			
Tungsten & compounds (as W) insoluble compounds & others			
	WEL short-term (15 min reference period)	10	mg/m ³
	WEL long-term (8-hr TWA reference period)	5	mg/m ³

8.2 Exposure controls**Appropriate engineering controls**

Ensure adequate ventilation, local exhaust at the work station if necessary.

Personal protective equipment**Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of dust formation, take appropriate measures for breathing protection in the event that workplace threshold values are not specified.

Eye / face protection

Safety glasses (EN 166)

Hand protection

Not necessary, if processed under regular conditions.

Other

Normal chemical work clothing.

Environmental exposure controls

Trade name: Tungsten

Current version : 4.0.2, issued: 03.04.2023

Replaced version: 4.0.1, issued: 27.10.2022

Region: GB

No data available.

SECTION 9: Physical and chemical properties**9.1 Information on basic physical and chemical properties**

State of aggregation			
solid			
Form			
solid			
Colour			
grey			
Odour			
odourless			
pH value			
No data available			
Boiling point / boiling range			
Value	appr.	5930	°C
Melting point/freezing point			
Value		3370	°C
Decomposition temperature			
No data available			
Flash point			
Not applicable			
Ignition temperature			
No data available			
Flammability			
No data available			
Lower explosion limit			
No data available			
Upper explosion limit			
No data available			
Vapour pressure			
No data available			
Relative vapour density			
No data available			
Relative density			
No data available			
Density			
No data available			
Bulk density			
Value	appr.	4000	kg/m³
Solubility in water			
Comments	insoluble		
Solubility			
No data available			
Partition coefficient n-octanol/water (log value)			
No data available			
Kinematic viscosity			

Trade name: Tungsten

Current version : 4.0.2, issued: 03.04.2023

Replaced version: 4.0.1, issued: 27.10.2022

Region: GB

No data available

Particle characteristics

No data available

9.2 Other information**Other information**

No data available.

SECTION 10: Stability and reactivity**10.1 Reactivity**

No data available.

10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

10.3 Possibility of hazardous reactions

Dangerous reactions are not to be expected when handling product according to its intended use.

10.4 Conditions to avoid

Extreme heat

10.5 Incompatible materials

metal oxides; Halogens; Halogenous-halogen compound; Fluorides; oxidizing agents; Peroxides

10.6 Hazardous decomposition products

No data available.

SECTION 11: Toxicological information**11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

No data available

Acute dermal toxicity

No data available

Acute inhalational toxicity

No data available

Skin corrosion/irritation

No data available

Serious eye damage/irritation

No data available

Respiratory or skin sensitisation

No data available

Germ cell mutagenicity

No data available

Reproduction toxicity

No data available

Carcinogenicity

No data available

STOT - single exposure

No data available

STOT - repeated exposure

No data available

Aspiration hazard

No data available

Trade name: Tungsten

Current version : 4.0.2, issued: 03.04.2023

Replaced version: 4.0.1, issued: 27.10.2022

Region: GB

Symptoms related to the physical, chemical and toxicological characteristics

Eye contact may cause mechanical irritation through dust particles.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Inhalation of high concentrations may result in damage of the lungs.

11.2 Information on other hazards**Endocrine disrupting properties**

No data available.

Other information

No data available.

SECTION 12: Ecological information**12.1 Toxicity****Toxicity to fish (acute)**

No data available

Toxicity to fish (chronic)

No data available

Toxicity to Daphnia (acute)

No data available

Toxicity to Daphnia (chronic)

No data available

Toxicity to algae (acute)

No data available

Toxicity to algae (chronic)

No data available

Bacteria toxicity

No data available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment**Results of PBT and vPvB assessment**

PBT assessment

The study does not need to be conducted according to Annex XIII of Regulation (EC) 1907/2006 (REACH).

vPvB assessment

The study does not need to be conducted according to Annex XIII of Regulation (EC) 1907/2006 (REACH).

12.6 Endocrine disrupting properties

No data available.

12.7 Other adverse effects

No data available.

12.8 Other information**Other information**

Do not discharge product unmonitored into the environment.

Ecological data are not available.

SECTION 13: Disposal considerations

Trade name: Tungsten**Current version :** 4.0.2, issued: 03.04.2023**Replaced version:** 4.0.1, issued: 27.10.2022**Region:** GB**13.1 Waste treatment methods****Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

Packaging

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse. Packaging that cannot be cleaned should be disposed of in agreement with the regional waste disposal company.

SECTION 14: Transport information**14.1 Transport ADR/RID/ADN**

The product is not subject to ADR/RID/ADN regulations.

14.2 Transport IMDG

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Maritime transport in bulk according to IMO instruments

Not relevant

SECTION 15: Regulatory information**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
EU regulations**Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

In accordance with the REACH regulation (EC) 1907/2006, the product does not contain any substances that are considered as subject to listing in annex XIV, inventory of substances requiring authorisation.

REACH candidate list of substances of very high concern (SVHC) for authorisation

In accordance with article 57 and article 59 of the Reach regulation (EC) 1907/2006, this substance is not considered as subject to listing in annex XIV, inventory of substances requiring authorisation ("Authorization list").

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES

The substance is not subject to the provisions of annex XVII (restriction entries) of the Reach regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This substance is not subject to Part 1 or 2 of Annex I

15.2 Chemical safety assessment

A chemical safety assessment has not been carried out for this substance.

SECTION 16: Other information**Further information**

Article no. 90220, 2.500 g

Sources of key data used to compile the data sheet:

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Trade name: Tungsten

Current version : 4.0.2, issued: 03.04.2023

Replaced version: 4.0.1, issued: 27.10.2022

Region: GB

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.
The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

Creation of the safety data sheet

UMCO GmbH

Georg-Wilhelm-Str. 187, D-21107 Hamburg

Tel.: +49 40 / 555 546 300 Fax: +49 40 / 555 546 357 e-mail: umco@umco.de

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Alterations/supplements:

Alterations to the previous edition are marked in the left-hand margin.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 640128